

INFO@GRUPOEMER.COM



INFO@GRUPOEMER.COM



Chemical	Letter	Level
40% Sodium hydroxide	К	2
Туре	С	

INTENT USE

The disposable nitrile gloves is intended to be worn on the hands of health care and similar personnel to prevent contamination between health care personnel and the patient's body. This is a single-use, powder-free, nonsterile device.

STANDARD COMPLIAINCE

PPE Cat III according to Regulation (EU) 2016/425

ISO 21420:2020 Protective gloves — General requirements and test methods EN ISO 374-1: 2016 Terminology and performance requirements for chemical risks EN 374-2:2014: Determination of resistance to penetration

EN 16523-1:2015+A1:2018 Permeation by potentially hazardous liquid chemicals under conditions of continuous contact EN ISO 374-4:2019 Determination of resistance to degradation by chemicals EN ISO 374-5:2016 Terminology and performance requirements for microorganisms risks

Medical Device Class I

EN 455-1: Requirements and testing for freedom from holes

EN 455-2: Requirements and testing for physical properties

EN 455-3: Requirements and testing for

biological evaluation

EN 455-4: Requirements and testing for shelf life determination

Food contact approved

FEATURE



Powder Free

Not made with natural rubber latex

Lab chemical tested

Dawn blue colour









